



## Product specification

### Cell culture flask, T-75, surface: Cell+, filter cap



## Product information

Product number	83.3911.302
----------------	-------------

## Product descriptions

Name	Cell culture flask, T-75, Surface: Cell+, Filter cap
Product number	83.3911.302
Product description	Cell culture flask, T-75, material: PS, surface: Cell+, for challenging adherent cells, coding colour: yellow, Filter cap, TC Tested, 5 piece(s)/bag

## Material & colours

Closure material	Polyethylene (PE)
Color of cap	yellow

## Product characteristics

Growth area	74.2 cm <sup>2</sup>
Vessel type	Flask
Cap	Filter cap
Cell type	challenging adherent cells

## Purity & certification

Product category	not regulated
CE certified	no CE mark

Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

## Packaging

Piece(s) / case	100
Piece(s) / pallet	1600
Height of case	438 mm
Width of case	402 mm
Depth of case	445 mm
Volume of case	0.0784 cbm
Weight of case	6.54 kg
EAN inner box	4038917337502
EAN case	4038917366076

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT