



Product specification

Cell culture flask, T-75, surface: Suspension, 2-position screw cap



Product information

Product number	83.3911.500
----------------	-------------

Product descriptions

Name	Cell culture flask, T-75, Surface: Suspension, 2-position screw cap
Product number	83.3911.500
Product description	Cell culture flask, T-75, material: PS, surface: Suspension, for suspension cells, coding colour: green, 2-position screw cap, TC Tested, 5 piece(s)/bag

Material & colours

Closure material	Polyethylene (PE)
Color of cap	green

Product characteristics

Growth area	74.2 cm ²
Vessel type	Flask
Cap	2-position screw cap
Cell type	suspension cells

Purity & certification

Product category	not regulated
------------------	---------------

CE certified	no CE mark
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / case	100
Piece(s) / pallet	1600
Height of case	438 mm
Width of case	402 mm
Depth of case	445 mm
Volume of case	0.0784 cbm
Weight of case	6.79 kg
EAN inner box	4038917337830
EAN case	4038917366083

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT