



## Product specification

### Cell culture flask, T-175, surface: Standard, 2-position screw cap



## Product information

Product number	83.3912
----------------	---------

## Product descriptions

Name	Cell culture flask, T-175, Surface: Standard, 2-position screw cap
Product number	83.3912
Product description	Cell culture flask, T-175, material: PS, surface: Standard, for adherent cells, coding colour: red, 2-position screw cap, TC Tested, 5 piece(s)/bag

## Material & colours

Product material	Polystyrene (PS)
Closure material	Polyethylene (PE)
Color of cap	red

## Product characteristics

Growth area	175 cm <sup>2</sup>
Vessel type	Flask
Cap	2-position screw cap
Cell type	adherent cells

## Purity & certification

Product category	no medical device   IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

## Packaging

Minimum order quantity	40
Piece(s) / case	40
Piece(s) / pallet	960
Height of case	450 mm
Width of case	267 mm
Depth of case	522 mm
Volume of case	0.0627 cbm
Weight of case	5.63 kg
EAN inner box	4038917337717
EAN case	4038917337700

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT