




## Product specification

### Cell culture plate, 6 well, surface: Standard, flat base



## Product information

Product number	83.3920.005
----------------	-------------

## Product descriptions

Name	Cell culture plate, 6 Well, Surface: Standard, Flat base
Product number	83.3920.005
Product description	Cell culture plate, 6 well, material: PS, surface: Standard, for adherent cells, coding colour: red, flat base, TC Tested, 5 piece(s)/bag

## Material & colours

Product material	Polystyrene (PS)
Color of product	transparent
Closure material	Polystyrene (PS)

## Size

Width of product	85.2 mm
Length of product	127.8 mm
Height of product	22.4 mm
Max. volume	15.75 ml

## Product characteristics

Growth area	9.07 cm <sup>2</sup>
Base shape	Flat base

Label/Print	With plastic label
Vessel type	Plate
Closure type	Slip-on lid
Cell type	adherent cells

## Purity & certification

Product category	no medical device   IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

## Packaging

Minimum order quantity	100
Piece(s) / case	100
Piece(s) / pallet	2000
Height of case	345 mm
Width of case	372 mm
Depth of case	567 mm
Volume of case	0.0728 cbm
Weight of case	7.26 kg
EAN inner box	4038917360135
EAN case	4038917360128

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT