


Product specification

Cell culture plate, 6 well, surface: Cell+, flat base



Product information

Order number	83.3920.300
--------------	-------------

Product descriptions

Name	Cell culture plate, 6 Well, Surface: Cell+, Flat base
Product number	83.3920.300
Product description	Cell culture plate, 6-well, PS, surface: Cell+ for sensitive adherent cells (e.g. primary cells), yellow, flat base, TC Tested (sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, DNase-/RNase-/DNA-free), 1 pc./paper/plastic packaging

Product characteristics

Growth area	8.87 cm ²
Base shape	Flat base
Label/Print	With plastic label
Tube	Plate
Type of cap	Slip-on lid
Cell type	Challenging, adherent

Size

Width of product	85.2 mm
Length of product	127.8 mm
Height of product	22.4 mm
Max. volume	15.75 ml

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

Cap material	Polystyrene (PS)
--------------	------------------

Purity & certification

Accordance with MDD or IVDD	no medical device IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Lot or serial number	yes

Packaging

Piece(s) / case	50
Piece(s) / pallet	2400
Height of case	275 mm
Width of case	297 mm
Depth of case	332 mm
Volume of case	0.0271 cbm
Weight of case	3 kg
EAN inner box	4038917338370
EAN case	4038917338363

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

