




## Product specification

### Cell culture plate, 12 well, surface: Cell+, flat base



## Product information

Product number	83.3921.300
----------------	-------------

## Product descriptions

Name	Cell culture plate, 12 Well, Surface: Cell+, Flat base
Product number	83.3921.300
Product description	Cell culture plate, 12 well, material: PS, surface: Cell+, for challenging adherent cells, coding colour: yellow, flat base, TC Tested, 1 piece(s)/blister

## Material & colours

Color of product	transparent
Closure material	Polystyrene (PS)

## Size

Width of product	85.2 mm
Length of product	127.8 mm
Height of product	22.5 mm
Max. volume	6.64 ml

## Product characteristics

Growth area	3.65 cm <sup>2</sup>
Number of wells	12
Base shape	Flat base

Label/Print	With plastic label
Vessel type	Plate
Closure type	Slip-on lid
Cell type	challenging adherent cells

## Purity & certification

Product category	no medical device   IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

## Packaging

Piece(s) / case	50
Piece(s) / pallet	2400
Height of case	275 mm
Width of case	297 mm
Depth of case	332 mm
Volume of case	0.0271 cbm
Weight of case	3.25 kg
EAN inner box	4038917338462
EAN case	4038917338455

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT