




Product specification

Cell culture plate, 24 well, surface: Standard, flat base



Product information

Product number	83.3922.005
----------------	-------------

Product descriptions

Name	Cell culture plate, 24 Well, Surface: Standard, Flat base
Product number	83.3922.005
Product description	Cell culture plate, 24 well, material: PS, surface: Standard, for adherent cells, coding colour: red, flat base, TC Tested, 5 piece(s)/bag

Material & colours

Color of product	transparent
Closure material	Polystyrene (PS)

Size

Width of product	85.2 mm
Length of product	127.8 mm
Height of product	22.45 mm
Max. volume	3.39 ml

Product characteristics

Growth area	1.82 cm ²
Number of wells	24
Base shape	Flat base

Label/Print	With plastic label
Vessel type	Plate
Closure type	Slip-on lid
Cell type	adherent cells

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / case	100
Piece(s) / pallet	2000
Height of case	335 mm
Width of case	372 mm
Depth of case	567 mm
Volume of case	0.0707 cbm
Weight of case	8.45 kg
EAN inner box	4038917360197
EAN case	4038917360180

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT