




Product specification

Cell culture plate, 96 well, surface: Suspension, flat base



Product information

Product number	83.3924.500
----------------	-------------

Product descriptions

Name	Cell culture plate, 96 Well, Surface: Suspension, Flat base
------	---

Product number	83.3924.500
----------------	-------------

Product description	Cell culture plate, 96 well, material: PS, surface: Suspension, for suspension cells, coding colour: green, flat base, TC Tested, 1 piece(s)/blister
---------------------	--

Material & colours

Product material	Polystyrene (PS)
------------------	------------------

Color of product	transparent
------------------	-------------

Closure material	Polystyrene (PS)
------------------	------------------

Size

Width of product	85.2 mm
------------------	---------

Length of product	127.8 mm
-------------------	----------

Height of product	16.55 mm
-------------------	----------

Max. volume	0.39 ml
-------------	---------

Product characteristics

Growth area	0.29 cm ²
-------------	----------------------

Number of wells	96
Base shape	Flat base
Label/Print	With plastic label
Vessel type	Plate
Closure type	Slip-on lid
Cell type	suspension cells

Purity & certification

Product category	no medical device IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	50
Piece(s) / case	50
Piece(s) / pallet	3200
Height of case	215 mm
Width of case	297 mm
Depth of case	332 mm
Volume of case	0.0212 cbm
Weight of case	3.54 kg
EAN inner box	4038917338769
EAN case	4038917338752

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT