



## Product specification

### Cell culture plate, 96 well, surface: Standard, round base



## Product information

Product number	83.3925
----------------	---------

## Product descriptions

Name	Cell culture plate, 96 Well, Surface: Standard, Round base
Product number	83.3925
Product description	Cell culture plate, 96 well, material: PS, surface: Standard, for adherent cells, coding colour: red, round base, TC Tested, 1 piece(s)/blister

## Material & colours

Color of product	transparent
Closure material	Polystyrene (PS)

## Size

Width of product	85.2 mm
Length of product	127.8 mm
Height of product	16.55 mm
Max. volume	0.31 ml

## Product characteristics

Number of wells	96
Base shape	Round base
Label/Print	With plastic label

Vessel type	Plate
Closure type	Slip-on lid
Cell type	adherent cells

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

## Packaging

Piece(s) / case	50
Piece(s) / pallet	3200
Height of case	205 mm
Width of case	297 mm
Depth of case	332 mm
Volume of case	0.0202 cbm
Weight of case	3.19 kg
EAN inner box	4038917338790
EAN case	4038917338783

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT