



## Product specification

**TC insert, for 6-well plates, PET, translucent, pore size: 5  $\mu\text{m}$**



## Product information

Product number	83.3930.500
----------------	-------------

## Product descriptions

Name	TC insert, for 6-well plates, PET, translucent, pore size: 5 $\mu\text{m}$
Product number	83.3930.500
Product description	TC insert, for 6-well plates, membrane: PET, translucent, pore size: 5 $\mu\text{m}$ , sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

## Material & colours

Membrane material	Polyethylene therephthalate (PET)
Frame material	Polystyrene (PS)

## Size

Height of product	16.38 mm
-------------------	----------

## Product characteristics

Growth area	4.5 cm <sup>2</sup>
Base shape	Flat base
Vessel type	TC insert
Pore size	5 $\mu\text{m}$

## Purity & certification

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	X-Ray
Pyrogen-free	yes
Batched	yes

## Packaging

Minimum order quantity	24
Piece(s) / inner box	24
Piece(s) / case	720
Piece(s) / pallet	11520
Height of inner box	131 mm
Width of inner box	104 mm
Depth of inner box	196 mm
Height of case	550 mm
Width of case	435 mm
Depth of case	410 mm
Volume of case	0.0981 cbm
Weight of case	6.4598 kg
EAN inner box	4038917370950
EAN case	4038917370943

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT