



Product specification

**TC insert, for 12-well plates, PET, transparent, pore size:
0.4 µm**



Product information

Product number	83.3931.041
----------------	-------------

Product descriptions

Name	TC insert, for 12-well plates, PET, transparent, pore size: 0.4 µm
Product number	83.3931.041
Product description	TC insert, for 12-well plates, membrane: PET, transparent, pore size: 0.4 µm, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Material & colours

Membrane material	Polyethylene therephthalate (PET)
Frame material	Polystyrene (PS)

Size

Height of product	16.38 mm
-------------------	----------

Product characteristics

Growth area	1.1 cm ²
Base shape	Flat base
Vessel type	TC insert

Pore size

0.4 µm

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	X-Ray
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	48
Piece(s) / case	1440
Piece(s) / pallet	23040
Height of inner box	131 mm
Width of inner box	104 mm
Depth of inner box	196 mm
Height of case	440 mm
Width of case	415 mm
Depth of case	550 mm
Volume of case	0.1004 cbm
Weight of case	8.5594 kg
EAN inner box	4038917371162
EAN case	4038917371155

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT