



## Product specification

**Filtropur V50, Vacuum filtration unit, 500 ml, PES, 0.45 µm**



### Product information

Product number	83.3941
----------------	---------

### Product descriptions

Name	Filtropur V50, Vacuum filtration unit, 500 ml, PES, 0.45 µm
Product number	83.3941
Product description	Filtropur V50, Vacuum filtration unit, 500 ml, screw cap, membrane: PES, Ø membraner: 90 mm, pore size: 0.45 µm, for clear filtration, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

### Material & colours

Color of product	transparent
Membrane material	Polyethersulfone (PES)
Closure material	Polypropylene (PP)

### Size

Height of product	21.84 cm
Sample volume	500 ml
Filtration volume	500 ml

### Product characteristics

Filter	yes
Area of the filter	63.6 cm <sup>2</sup>

Graduation	yes
Pore size	0.45 µm
Closure type	Screw cap
Connection	Step taper
Example of use	for clear filtration

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

## Packaging

Piece(s) / case	12
Piece(s) / pallet	288
Height of inner box	220 mm
Width of inner box	115 mm
Depth of inner box	115 mm
Height of case	265 mm
Width of case	375 mm
Depth of case	527 mm
Volume of case	0.0524 cbm
Weight of case	3.36 kg
EAN inner box	4038917412025
EAN case	4038917412018

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT