



## Product specification

### Filtropur V50, Vacuum filtration unit, 500 ml, PES, 0.2 µm



## Product information

Product number	83.3941.001
----------------	-------------

## Product descriptions

Name	Filtropur V50, Vacuum filtration unit, 500 ml, PES, 0.2 µm
Product number	83.3941.001
Product description	Filtropur V50, Vacuum filtration unit, 500 ml, screw cap, membrane: PES, Ø membraner: 90 mm, pore size: 0.2 µm, for sterile filtration, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

## Material & colours

Product material	Polystyrene (PS)
Color of product	transparent
Membrane material	Polyethersulfone (PES)
Closure material	Polypropylene (PP)

## Size

Height of product	21.84 cm
Sample volume	500 ml
Filtration volume	500 ml

## Product characteristics

Filter	yes
--------	-----

Area of the filter	63.6 cm <sup>2</sup>
Graduation	yes
Pore size	0.2 µm
Closure type	Screw cap
Connection	Step taper
Example of use	for sterile filtration

## Purity & certification

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

## Packaging

Minimum order quantity	12
Piece(s) / case	12
Piece(s) / pallet	288
Height of case	381 mm
Width of case	267 mm
Depth of case	508 mm
Volume of case	0.0516 cbm
Weight of case	3.36 kg
EAN inner box	4038917412056
EAN case	4038917412049

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT