



## Product specification

### Filtropur BT 50, Bottle top filter, 500 ml, PES, 0.45 µm



## Product information

Product number	83.3941.100
----------------	-------------

## Product descriptions

Name	Filtropur BT 50, Bottle top filter, 500 ml, PES, 0.45 µm
Product number	83.3941.100
Product description	Filtropur BT 50, Bottle top filter, for screw cap 45 mm, 500 ml, membrane: PES, Ø membraner: 90 mm, pore size: 0.45 µm, for clear filtration, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

## Material & colours

Product material	Polystyrene (PS)
Color of product	transparent
Membrane material	Polyethersulfone (PES)

## Size

Height of product	10.61 cm
Sample volume	500 ml
Filtration volume	500 ml

## Product characteristics

Filter	yes
Area of the filter	63.6 cm <sup>2</sup>

Graduation	yes
Pore size	0.45 µm
Connection	Step taper
Example of use	for clear filtration

## Purity & certification

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

## Packaging

Minimum order quantity	24
Piece(s) / case	24
Piece(s) / pallet	576
Height of case	381 mm
Width of case	267 mm
Depth of case	521 mm
Volume of case	0.0529 cbm
Weight of case	3.49 kg
EAN inner box	4038917412230
EAN case	4038917412223

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT