




Product specification

Cell strainer, pore size: 40 µm, blue



Product information

| | |
|----------------|-------------|
| Product number | 83.3945.040 |
|----------------|-------------|

Product descriptions

| | |
|---------------------|---|
| Name | Cell strainer, Pore size: 40 µm, blue |
| Product number | 83.3945.040 |
| Product description | Cell strainer, pore size: 40 µm, blue, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister |

Material & colours

| | |
|-------------------|--------------------|
| Color of product | blue |
| Filter material | Nylon |
| Membrane material | Nylon |
| Frame material | Polypropylene (PP) |
| Color of frame | blue |

Size

| | |
|-------------------|----------|
| Diameter | 34.09 mm |
| Height of product | 21 mm |

Product characteristics

| | |
|--------------------|---------------------|
| Filter | yes |
| Area of the filter | 4.2 cm ² |
| Pore size | 40 µm |

Purity & certification

| | |
|------------------|--|
| Product category | no medical device IVD |
| Purity standard | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation | Electron irradiation |
| Pyrogen-free | yes |
| Batched | yes |

Packaging

| | |
|------------------------|---------------|
| Minimum order quantity | 50 |
| Piece(s) / inner box | 50 |
| Piece(s) / case | 600 |
| Piece(s) / pallet | 9600 |
| Height of case | 410 mm |
| Width of case | 385 mm |
| Depth of case | 520 mm |
| Volume of case | 0.0821 cbm |
| Weight of case | 6.78 kg |
| EAN inner box | 4038917414982 |
| EAN case | 4038917414975 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT