




## Product specification

**Cell strainer, pore size: 40 µm, blue**



### Product information

Product number	83.3945.040
----------------	-------------

### Product descriptions

Name	Cell strainer, Pore size: 40 µm, blue
Product number	83.3945.040
Product description	Cell strainer, pore size: 40 µm, blue, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

### Material & colours

Color of product	blue
Filter material	Nylon
Membrane material	Nylon
Frame material	Polypropylene (PP)
Color of frame	blue

### Size

Diameter	34.09 mm
Height of product	21 mm

### Product characteristics

Filter	yes
Area of the filter	4.2 cm <sup>2</sup>
Pore size	40 µm

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

## Packaging

Piece(s) / inner box	50
Piece(s) / case	600
Piece(s) / pallet	9600
Height of inner box	130 mm
Width of inner box	180 mm
Depth of inner box	245 mm
Height of case	410 mm
Width of case	385 mm
Depth of case	520 mm
Volume of case	0.0821 cbm
Weight of case	6.78 kg
EAN inner box	4038917414982
EAN case	4038917414975

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT