




Product specification

Cell scraper, 2-position blade, size: S



Product information

Product number	83.3950
----------------	---------

Product descriptions

Name	Cell scraper, 2-position blade, size: S
Product number	83.3950
Product description	Cell scraper, 2-position blade, size: S, handle length: 24 cm, blade length: 1.35 cm, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Material & colours

Product material	Polystyrene (PS)
Color of product	red/opaque
Blade material	Thermoplastic elastomers (TPE)
Stem material	Polystyrene (PS)

Size

Lennght of blade	1.35 cm
Length of product	24.7 cm
Width of blade	1.35 cm
Length of the stem	24 cm

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	100
Piece(s) / case	100
Piece(s) / pallet	16500
Height of case	333 mm
Width of case	110 mm
Depth of case	276 mm
Volume of case	0.0101 cbm
Weight of case	1.3 kg
EAN inner box	4038917426633
EAN case	4038917426626

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT