




Product specification

Cell scraper, 2-position blade, size: M



Product information

Product number	83.3951
----------------	---------

Product descriptions

Name	Cell scraper, 2-position blade, size: M
Product number	83.3951
Product description	Cell scraper, 2-position blade, size: M, handle length: 24 cm, blade length: 1.7 cm, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Material & colours

Color of product	red/opaque
Blade material	Thermoplastic elastomers (TPE)
Stem material	Polystyrene (PS)

Size

Lennght of blade	1.7 cm
Length of product	24.9 cm
Width of blade	1.7 cm
Length of the stem	24 cm

Purity & certification

Product category	not regulated
------------------	---------------

CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / case	100
Piece(s) / pallet	16500
Height of inner box	7 mm
Width of inner box	20 mm
Depth of inner box	247 mm
Height of case	110 mm
Width of case	270 mm
Depth of case	330 mm
Volume of case	0.0098 cbm
Weight of case	1.3 kg
EAN inner box	4038917426664
EAN case	4038917426657

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT