




## Product specification

### Cell scraper, 2-position blade, size: L



## Product information

Product number	83.3952
----------------	---------

## Product descriptions

Name	Cell scraper, 2-position blade, size: L
Product number	83.3952
Product description	Cell scraper, 2-position blade, size: L, handle length: 36 cm, blade length: 1.7 cm, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

## Material & colours

Color of product	red/opaque
Blade material	Thermoplastic elastomers (TPE)
Stem material	Polystyrene (PS)

## Size

Lennght of blade	1.7 cm
Length of product	36.9 cm
Width of blade	1.7 cm
Length of the stem	36 cm

## Purity & certification

Product category	not regulated
------------------	---------------

CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

## Packaging

Piece(s) / case	100
Piece(s) / pallet	7000
Height of inner box	10 mm
Width of inner box	20 mm
Depth of inner box	365 mm
Height of case	137 mm
Width of case	340 mm
Depth of case	495 mm
Volume of case	0.0231 cbm
Weight of case	2.2 kg
EAN inner box	4038917426695
EAN case	4038917426688

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT