



Product specification

Safety lancet, Mini, Ø needle: 28 G, penetration depth: 1.6 mm



Product information

Product number	85.1015
----------------	---------

Product descriptions

Name	Safety lancet, Mini, Ø needle: 28 G, Penetration depth: 1.6 mm
Product number	85.1015
Product description	Safety lancet, Mini, Ø needle: 28 G, light blue, penetration depth: 1.6 mm, for capillary blood collection, for puncturing finger, ear, heel, 200 piece(s)/case

Material & colours

Color of product	light blue
------------------	------------

Size

Width of product	21 mm
Length of product	52 mm
Diameter of needle (external)	28 G
Blood volume	0.1 µl - 10 µl

Product characteristics

Safety product	yes
----------------	-----

Example of use

Young children, adults: diabetic patients, oncology, determining blood sugar

Penetration depth

1.6 mm

Puncture site

Finger, ear, heel

Purity & certification

Product category

Medical device (MD)

CE certified

CE0124 - DEKRA

Purity standard

sterile

Sterilisation

Gamma irradiation

Batched

yes

Packaging

Piece(s) / inner box

200

Piece(s) / case

1200

Piece(s) / pallet

108000

Height of inner box

151 mm

Width of inner box

118 mm

Depth of inner box

118 mm

Height of case

173 mm

Width of case

248 mm

Depth of case

366 mm

Volume of case

0.0157 cbm

Weight of case

4.08 kg

EAN inner box

4038917085021

EAN case

4038917057622

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT