



Product specification

**Safety lancet, Extra, Ø needle: 18 G, penetration depth:
1.8 mm**



Product information

Product number	85.1017
----------------	---------

Product descriptions

Name	Safety lancet, Extra, Ø needle: 18 G, Penetration depth: 1.8 mm
Product number	85.1017
Product description	Safety lancet, Extra, Ø needle: 18 G, yellow, penetration depth: 1.8 mm, for capillary blood collection, for puncturing finger, ear, heel, 200 piece(s)/case

Material & colours

Color of product	yellow
------------------	--------

Size

Width of product	21 mm
Length of product	52 mm
Diameter of needle (external)	18 G
Blood volume	20 µl - 100 µl

Product characteristics

Safety product	yes
----------------	-----

Example of use	Children, Adults, obese adults: capillary blood gas analysis, CBC
Penetration depth	1.8 mm
Puncture site	Finger, ear, heel

Purity & certification

Product category	Medical device
CE certified	CE 0124 (DEKRA)
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order quantity	1200
Piece(s) / inner box	200
Piece(s) / case	1200
Piece(s) / pallet	108000
Height of inner box	151 mm
Width of inner box	118 mm
Depth of inner box	118 mm
Height of case	173 mm
Width of case	248 mm
Depth of case	366 mm
Volume of case	0.0157 cbm
Weight of case	4.08 kg
EAN inner box	4038917085045
EAN case	4038917057646

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT