



## Product specification

**S-Monovette® needle, 20G x 1 1/2", yellow, 1 piece(s)/blister**



## Product information

Product number	85.1160
----------------	---------

## Product descriptions

Name	S-Monovette® needle, 20G x 1 1/2", yellow, 1 piece(s)/blister
Product number	85.1160
Product description	S-Monovette® needle, needle (ØxL): 20G x 1 1/2", coding colour: yellow, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 100 piece(s)/case

## Product characteristics

Safety product	no
Connection	Membrane screw cap

## Purity & certification

Product category	Medical device (MD)
CE certified	CE0124 - DEKRA
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gamma irradiation
Pyrogen-free	yes

Batched	yes
Veterinary Use Only	yes

## Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	78000
Height of inner box	71 mm
Width of inner box	111 mm
Depth of inner box	268 mm
Height of case	290 mm
Width of case	233 mm
Depth of case	373 mm
Volume of case	0.0252 cbm
Weight of case	2.96 kg
EAN inner box	4038917085519
EAN case	4038917004169

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT