



Product specification

S-Monovette® safety needle, 21G x 1 1/2", green, 1 piece(s)/blister



Product information

| | |
|----------------|-------------|
| Product number | 85.1162.600 |
|----------------|-------------|

Product descriptions

| | |
|---------------------|--|
| Name | S-Monovette® safety needle, 21G x 1 1/2", green, 1 piece(s)/blister |
| Product number | 85.1162.600 |
| Product description | S-Monovette® safety needle, for S-Monovette® metal analysis, needle (ØxL): 21G x 1 1/2", coding colour: green, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 50 piece(s)/case |

Product characteristics

| | |
|----------------|--------------------|
| Safety product | yes |
| Connection | Membrane screw cap |

Purity & certification

| | |
|------------------|---------------------------------------|
| Product category | Medical device (MD) |
| CE certified | CE0124 - DEKRA |
| Purity standard | sterile, non-pyrogenic/endotoxin-free |
| Sterilisation | Gamma irradiation |
| Pyrogen-free | yes |

Batched

yes

Packaging

| | |
|----------------------|---------------|
| Piece(s) / inner box | 50 |
| Piece(s) / case | 500 |
| Piece(s) / pallet | 39000 |
| Height of inner box | 71 mm |
| Width of inner box | 111 mm |
| Depth of inner box | 268 mm |
| Height of case | 290 mm |
| Width of case | 233 mm |
| Depth of case | 373 mm |
| Volume of case | 0.0252 cbm |
| Weight of case | 2.48 kg |
| EAN inner box | 4038917250474 |
| EAN case | 4038917361934 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT