


Product specification

Multifly® needle, 21G x 3/4", green, tube length: 80 mm, 1 piece(s)/blister



Product information

Product number	85.1638.005
----------------	-------------

Product descriptions

Name	Multifly® needle, 21G x 3/4", green, Tube length: 80 mm, 1 piece(s)/blister
------	--

Product number	85.1638.005
----------------	-------------

Product description	Multifly® needle, needle (ØxL): 21G x 3/4", coding colour: green, tube length: 80 mm, with Multi-adapter, DEHP-free, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 120 piece(s)/case
---------------------	---

Size

Diameter of tube	2.3 mm
------------------	--------

Length of tube	80 mm
----------------	-------

Product characteristics

Design	Multi-adapter
--------	---------------

Safety product	no
----------------	----

Connection	Membrane screw cap
------------	--------------------

Product category	Medical device (MD)
CE certified	CE0124 - DEKRA
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Plasticiser	DEHP-free
Batched	yes
Veterinary Use Only	yes

Packaging

Piece(s) / inner box	120
Piece(s) / case	480
Piece(s) / pallet	20640
Height of inner box	171 mm
Width of inner box	214 mm
Depth of inner box	224 mm
Height of case	240 mm
Width of case	358 mm
Depth of case	442 mm
Volume of case	0.038 cbm
Weight of case	3.1901 kg
EAN inner box	4038917423762
EAN case	4038917423472

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

