



## Product specification

### miniPERM®, Start-up support kit, for miniPERM® bioreactor



## Product information

Product number	94.6001.094
----------------	-------------

## Product descriptions

Name	miniPERM®, Start-up support kit, for miniPERM® bioreactor
Product number	94.6001.094
Product description	miniPERM®, Start-up support kit, accessories for inoculation, sample collection and harvesting, for miniPERM® bioreactor, not autoclavable, set components: 8 x disposable syringe 50 ml, 20 x disposable syringe 2 ml, 8 x filling tube, 20 x needle, 6 x septum cap, 1 x miniPERM® stand, 1 piece(s)/bag

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

## Packaging

Piece(s) / pallet	96
-------------------	----

Height of inner box	80 mm
Width of inner box	150 mm
Depth of inner box	200 mm
Height of case	140 mm
Width of case	278 mm
Depth of case	348 mm
Volume of case	0.0135 cbm
Weight of case	0.93 kg
EAN inner box	4038917254939
EAN case	4038917254922

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT