




Product specification

**Safety-Multifly® needle, 21G x 3/4", green, tube length:
200 mm, 1 piece(s)/blister**



Product information

Product number	85.1638.200
----------------	-------------

Product descriptions

Name	Safety-Multifly® needle, 21G x 3/4", green, Tube length: 200 mm, 1 piece(s)/blister
Product number	85.1638.200
Product description	Safety-Multifly® needle, needle (ØxL): 21G x 3/4", coding colour: green, tube length: 200 mm, with Multi-Adapter, screw-on guide sleeve, DEHP-free, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 120 piece(s)/case

Size

Diameter of tube	2.3 mm
Length of tube	200 mm

Product characteristics

Design	Multi-Adapter, screw-on guide sleeve
Safety product	yes
Connection	Membrane screw cap

Purity & certification

Product category	Medical device (MD)
CE certified	CE0124 - DEKRA
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Plasticiser	DEHP-free
Batched	yes

Packaging

Piece(s) / inner box	120
Piece(s) / case	480
Piece(s) / pallet	20640
Height of inner box	171 mm
Width of inner box	214 mm
Depth of inner box	224 mm
Height of case	240 mm
Width of case	358 mm
Depth of case	442 mm
Volume of case	0.038 cbm
Weight of case	3.66 kg
EAN inner box	4038917423885
EAN case	4038917423595

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT