




Product specification

Micro needle, 21G x 3/4", green, 1 piece(s)/blister



Product information

Product number	85.921
----------------	--------

Product descriptions

Name	Micro needle, 21G x 3/4", green, 1 piece(s)/blister
Product number	85.921
Product description	Micro needle, for collecting very small amounts of blood, needle (ØxL): 21G x 3/4", coding colour: green, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 50 piece(s)/case

Product characteristics

Safety product	no
----------------	----

Purity & certification

Product category	Medical device
CE certified	CE 0124 (DEKRA)
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	50
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	55000
Height of inner box	66 mm
Width of inner box	108 mm
Depth of inner box	139 mm
Height of case	160 mm
Width of case	228 mm
Depth of case	338 mm
Volume of case	0.0123 cbm
Weight of case	1.37 kg
EAN inner box	4038917315517
EAN case	4038917315500

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT