



## Product specification

# Universal turning device, for 4 miniPERM® bioreactors, 115/230 V



## Product information

Product number	94.6001.061
----------------	-------------

## Product descriptions

Name	Universal turning device, for 4 miniPERM® bioreactors, 115/230 V
Product number	94.6001.061
Product description	Universal turning device, for 4 miniPERM® bioreactors, not autoclavable, 115/230 V, 1 piece(s)/case

## Size

Width of product	26 cm
Length of product	46 cm

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

## Packaging

Piece(s) / pallet	30
-------------------	----

Height of inner box	325 mm
Width of inner box	358 mm
Depth of inner box	398 mm
Height of case	325 mm
Width of case	358 mm
Depth of case	398 mm
Volume of case	0.0463 cbm
Weight of case	6.3 kg
EAN inner box	4038917254878
EAN case	4038917254861

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT