




## Product specification

**Transfer pipette, 3.5 ml, (LxW): 156 x 12.5 mm, LD-PE,  
transparent**



## Product information

Product number	86.1172.020
----------------	-------------

## Product descriptions

Name	Transfer pipette, 3.5 ml, (LxW): 156 x 12.5 mm, LD-PE, transparent
Product number	86.1172.020
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 156 x 12.5 mm, material: LD-PE, transparent, droplet size: 30 µl - 45 µl, sterile, 20 piece(s)/bag

## Material & colours

Product material	Low Density Polyethylene (LD-PE)
Color of product	transparent

## Size

Width of product	12.5 mm
Length of product	156 mm
Sample volume	3.5 ml
Droplet size min-max	30 µl - 45 µl

## Product characteristics

Graduation	yes
------------	-----

## Purity & certification

Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Minimum order quantity	4000
Piece(s) / inner box	500
Piece(s) / case	4000
Piece(s) / pallet	64000
Volume of case	0.1309 cbm
Weight of case	8.64 kg
EAN inner box	4038917318310
EAN case	4038917318303

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT