




Product specification

Transfer pipette, 3.5 ml, (LxW): 156 x 12.5 mm, LD-PE, transparent



Product information

Product number	86.1172.301
----------------	-------------

Product descriptions

Name	Transfer pipette, 3.5 ml, (LxW): 156 x 12.5 mm, LD-PE, transparent
Product number	86.1172.301
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 156 x 12.5 mm, material: LD-PE, transparent, droplet size: 30 µl - 45 µl, sterile, 1 piece(s)/blister

Material & colours

Color of product	transparent
------------------	-------------

Size

Width of product	12.5 mm
Length of product	156 mm
Sample volume	3.5 ml
Droplet size min-max	30 µl - 45 µl

Product characteristics

Graduation	yes
------------	-----

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	4000
Piece(s) / pallet	30240
Volume of case	0.1359 cbm
Weight of case	13.092 kg
EAN inner box	4038917378444
EAN case	4038917378437

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT