



Product specification

**Transfer pipette, 6 ml, (LxW): 146 x 15 mm, LD-PE,
transparent**



Product information

Product number	86.1175.020
----------------	-------------

Product descriptions

Name	Transfer pipette, 6 ml, (LxW): 146 x 15 mm, LD-PE, transparent
Product number	86.1175.020
Product description	Transfer pipette, nominal volume: 6 ml, (LxW): 146 x 15 mm, material: LD-PE, transparent, droplet size: 20 µl - 30 µl, sterile, 20 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Width of product	15 mm
Length of product	146 mm
Sample volume	6 ml
Droplet size min-max	20 µl - 30 µl

Product characteristics

Graduation	no
------------	----

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Piece(s) / inner box	400
Piece(s) / case	3200
Piece(s) / pallet	3200
Volume of case	0.1309 cbm
Weight of case	8.64 kg
EAN inner box	4038917318433
EAN case	4038917318426

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT