




## Product specification

**Transfer pipette, 6 ml, (LxW): 146 x 15 mm, LD-PE,  
transparent**



## Product information

Product number	86.1175.300
----------------	-------------

## Product descriptions

Name	Transfer pipette, 6 ml, (LxW): 146 x 15 mm, LD-PE, transparent
Product number	86.1175.300
Product description	Transfer pipette, nominal volume: 6 ml, (LxW): 146 x 15 mm, material: LD-PE, transparent, droplet size: 20 µl - 30 µl, 400 piece(s)/case

## Material & colours

Color of product	transparent
------------------	-------------

## Size

Width of product	15 mm
Length of product	146 mm
Sample volume	6 ml
Droplet size min-max	20 µl - 30 µl

## Product characteristics

Graduation	no
------------	----

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

## Packaging

Piece(s) / inner box	400
Piece(s) / case	4000
Piece(s) / pallet	44800
Volume of case	0.0994 cbm
Weight of case	7.952 kg
EAN inner box	4038917378536
EAN case	4038917378529

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT