




Product specification

**Transfer pipette, 2 ml, (LxW): 154 x 11 mm, LD-PE,
transparent**



Product information

Product number	86.1176.300
----------------	-------------

Product descriptions

Name	Transfer pipette, 2 ml, (LxW): 154 x 11 mm, LD-PE, transparent
Product number	86.1176.300
Product description	Transfer pipette, nominal volume: 2 ml, (LxW): 154 x 11 mm, material: LD-PE, transparent, droplet size: 35 µl - 45 µl

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Color of product	transparent

Size

Width of product	11 mm
Length of product	154 mm
Sample volume	2 ml
Droplet size min-max	35 µl - 45 µl

Product characteristics

Graduation	no
------------	----

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order quantity	5000
Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	84000
Volume of case	0.0994 cbm
Weight of case	8.64 kg
EAN inner box	4038917378567
EAN case	4038917378550

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT