



## Product specification

**Disposal bags, 24 l, (LxW): 780 x 400 mm, PP, transparent, without print**



### Product information

|                |         |
|----------------|---------|
| Product number | 86.1199 |
|----------------|---------|

### Product descriptions

|                     |   |
|---------------------|---|
| Name                | Disposal bags, 24 l, (LxW): 780 x 400 mm, PP, transparent, without print  |
| Product number      | 86.1199   |
| Product description | Disposal bags, volume: 24 l, (LxW): 780 x 400 mm, material: PP, transparent, film thickness: 50 µm, autoclavable up to 134 °C, without print, 50 piece(s)/bag |

### Material & colours

|                  |                    |
|------------------|--------------------|
| Product material | Polypropylene (PP) |
| Color of product | transparent        |

### Size

|                   |        |
|-------------------|--------|
| Width of product  | 400 mm |
| Length of product | 780 mm |
| Foil thickness    | 50 µm  |

### Product characteristics

|             |               |
|-------------|---------------|
| Label/Print | without print |
|-------------|---------------|

|                        |        |
|------------------------|--------|
| Autoclavable (max. °C) | 134 °C |
| Liquid-tight           | yes    |

## Purity & certification

|                   |                         |
|-------------------|-------------------------|
| Product category  | no medical device   IVD |
| Latex-free        | yes                     |
| Batched           | no                      |
| Safety data sheet | no                      |

## Packaging

|                        |               |
|------------------------|---------------|
| Minimum order quantity | 250           |
| Piece(s) / inner box   | 50            |
| Piece(s) / case        | 250           |
| Piece(s) / pallet      | 16500         |
| Height of case         | 155 mm        |
| Width of case          | 268 mm        |
| Depth of case          | 438 mm        |
| Volume of case         | 0.0182 cbm    |
| Weight of case         | 7.05 kg       |
| EAN inner box          | 4038917085731 |
| EAN case               | 4038917004367 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT