



Product specification

Disposal bags, 7 l, (LxW): 500 x 300 mm, PP, transparent, with print Biohazard



Product information

Product number	86.1201
----------------	---------

Product descriptions

Name	Disposal bags, 7 l, (LxW): 500 x 300 mm, PP, transparent, with print Biohazard
Product number	86.1201
Product description	Disposal bags, volume: 7 l, (LxW): 500 x 300 mm, material: PP, transparent, film thickness: 50 µm, autoclavable up to 134 °C, with print, print Biohazard, label/print: red, 50 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Width of product	300 mm
Length of product	500 mm
Foil thickness	50 µm

Product characteristics

Label/Print	With print
Color of print/label	red

Autoclavable (max. °C)	134 °C
Liquid-tight	yes

Purity & certification

Product category	not regulated
CE certified	no CE mark
Latex-free	yes
Batched	no
Safety data sheet	no

Packaging

Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	25000
Height of inner box	25 mm
Width of inner box	135 mm
Depth of inner box	305 mm
Height of case	140 mm
Width of case	278 mm
Depth of case	348 mm
Volume of case	0.0135 cbm
Weight of case	7.3 kg
EAN inner box	4038917085755
EAN case	4038917004381

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT