




## Product specification

**Disposal bags, 24 l, (LxW): 780 x 400 mm, PP, transparent, with print Biohazard**



## Product information

Product number	86.1202
----------------	---------

## Product descriptions

Name	Disposal bags, 24 l, (LxW): 780 x 400 mm, PP, transparent, with print Biohazard
Product number	86.1202
Product description	Disposal bags, volume: 24 l, (LxW): 780 x 400 mm, material: PP, transparent, film thickness: 50 µm, autoclavable up to 134 °C, with print, print Biohazard, label/print: red, 50 piece(s)/bag

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

## Size

Width of product	400 mm
Length of product	780 mm
Foil thickness	50 µm

## Product characteristics

Label/Print	With print
Color of print/label	red
Autoclavable (max. °C)	134 °C
Liquid-tight	yes

## Purity & certification

Product category	no medical device   IVD
Latex-free	yes
Batched	no
Safety data sheet	no

## Packaging

Minimum order quantity	250
Piece(s) / inner box	50
Piece(s) / case	250
Piece(s) / pallet	9000
Height of case	155 mm
Width of case	268 mm
Depth of case	438 mm
Volume of case	0.0182 cbm
Weight of case	7.1 kg
EAN inner box	4038917085762
EAN case	4038917004398

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT