




Product specification

**Serological pipette, with tip, unplugged, 1 ml, non-sterile,
200 piece(s)/bag**



Product information

Product number	86.1251
----------------	---------

Product descriptions

Name	Serological pipette, with tip, unplugged, 1 ml, non-sterile, 200 piece(s)/bag
Product number	86.1251
Product description	Serological pipette, with tip, unplugged, 1 ml, coding colour: yellow, material: PS, non-sterile, 200 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	4 mm
Length of product	270 mm
Sample volume	1 ml

Product characteristics

Design	with tip, unplugged
Label/Print	with printed graduations

Color of print/label	yellow/black
Graduation	yes

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	non-sterile
Batched	yes

Packaging

Piece(s) / inner box	200
Piece(s) / case	1000
Piece(s) / pallet	99000
Height of inner box	45 mm
Width of inner box	110 mm
Depth of inner box	290 mm
Height of case	140 mm
Width of case	250 mm
Depth of case	410 mm
Volume of case	0.0132 cbm
Weight of case	3.75 kg
EAN inner box	4038917085878
EAN case	4038917004497

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT