




## Product specification

**Serological pipette, with tip, unplugged, 1 ml, 200 piece(s)/bag**



### Product information

Product number	86.1251
----------------	---------

### Product descriptions

Name	Serological pipette, with tip, unplugged, 1 ml, 200 piece(s)/bag
Product number	86.1251
Product description	Serological pipette, with tip, unplugged, 1 ml, coding colour: yellow, material: PS, 200 piece(s)/bag

### Material & colours

Product material	Polystyrene (PS)
------------------	------------------

### Size

Diameter	4 mm
Length of product	270 mm

### Product characteristics

Design	with tip, unplugged
Color of print/label	yellow
Graduation	yes

### Purity & certification

Product category	no medical device   IVD
Batched	yes

## Packaging

Minimum order quantity	1000
Piece(s) / inner box	200
Piece(s) / case	1000
Piece(s) / pallet	99000
Volume of case	0.0132 cbm
Weight of case	3.75 kg
EAN inner box	4038917085878
EAN case	4038917004497

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT