




## Product specification

**Serological pipette, with tip, plugged, 2 ml, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic,**



## Product information

|                |             |
|----------------|-------------|
| Product number | 86.1252.001 |
|----------------|-------------|

## Product descriptions

|                     |  |
|---------------------|--|
| Name                | Serological pipette, with tip, plugged, 2 ml, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister                                     |
| Product number      | 86.1252.001  |
| Product description | Serological pipette, with tip, plugged, 2 ml, coding colour: green, material: PS, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister |

## Material & colours

|                  |             |
|------------------|-------------|
| Color of product | transparent |
|------------------|-------------|

## Size

|                   |        |
|-------------------|--------|
| Diameter          | 6 mm   |
| Length of product | 270 mm |
| Sample volume     | 2 ml   |

## Product characteristics

|             |                          |
|-------------|--------------------------|
| Design      | with tip, plugged        |
| Label/Print | with printed graduations |

|                      |             |
|----------------------|-------------|
| Color of print/label | green/black |
| Graduation           | yes         |

## Purity & certification

|                  |  |
|------------------|--|
| Product category | not regulated  |
| CE certified     | no CE mark   |
| Purity standard  | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation    | Gamma irradiation                                    |
| Pyrogen-free     | yes  |
| Batched          | yes  |

## Packaging

|                      |               |
|----------------------|---------------|
| Piece(s) / inner box | 100           |
| Piece(s) / case      | 1000          |
| Piece(s) / pallet    | 54000         |
| Height of case       | 290 mm        |
| Width of case        | 250 mm        |
| Depth of case        | 420 mm        |
| Volume of case       | 0.0323 cbm    |
| Weight of case       | 5.03 kg       |
| EAN inner box        | 4038917127196 |
| EAN case             | 4038917021081 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT