




Product specification

Serological pipette, with tip, unplugged, 5 ml, 100 piece(s)/bag



Product information

Product number	86.1253
----------------	---------

Product descriptions

Name	Serological pipette, with tip, unplugged, 5 ml, 100 piece(s)/bag
Product number	86.1253
Product description	Serological pipette, with tip, unplugged, 5 ml, coding colour: blue, material: PS, 100 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

Size

Diameter	8 mm
Length of product	340 mm

Product characteristics

Design	with tip, unplugged
Color of print/label	blue
Graduation	yes

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order quantity	500
Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	40500
Volume of case	0.0194 cbm
Weight of case	3.71 kg
EAN inner box	4038917085892
EAN case	4038917004510

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT