




Product specification

**Serological pipette, with tip, plugged, 10 ml, non-sterile,
100 piece(s)/bag**



Product information

Product number	86.1254
----------------	---------

Product descriptions

Name	Serological pipette, with tip, plugged, 10 ml, non-sterile, 100 piece(s)/bag
Product number	86.1254
Product description	Serological pipette, with tip, plugged, 10 ml, coding colour: orange, material: PS, non-sterile, 100 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	10 mm
Length of product	335 mm
Sample volume	10 ml

Product characteristics

Design	with tip, plugged
Label/Print	with printed graduations
Color of print/label	orange/black
Graduation	yes

Purity & certification

Product category	no medical device IVD
Purity standard	non-sterile
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	27000
Height of inner box	60 mm
Width of inner box	165 mm
Depth of inner box	340 mm
Height of case	260 mm
Width of case	250 mm
Depth of case	410 mm
Volume of case	0.0264 cbm
Weight of case	5.04 kg
EAN inner box	4038917085908
EAN case	4038917004527

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT