



Product specification

miniPERM®, Rack, for miniPERM® bioreactor



Product information

Product number	94.6001.054
----------------	-------------

Product descriptions

Name	miniPERM®, Rack, for miniPERM® bioreactor
Product number	94.6001.054
Product description	miniPERM®, Rack, for miniPERM® bioreactor, not autoclavable, 4 piece(s)/case

Material & colours

Product material	Acrylnitrile butadiene styrene (ABS)
Color of product	white

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Piece(s) / case	24
Piece(s) / pallet	672
Height of case	255 mm
Width of case	398 mm
Depth of case	588 mm

Volume of case	0.0597 cbm
Weight of case	3.7 kg
EAN inner box	4038917254786
EAN case	4038917254779

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT