




Product specification

Serological pipette, with tip, plugged, 50 ml, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic,



Product information

Product number	86.1256.001
----------------	-------------

Product descriptions

Name	Serological pipette, with tip, plugged, 50 ml, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister
Product number	86.1256.001
Product description	Serological pipette, with tip, plugged, 50 ml, coding colour: violet, material: PS, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	18 mm
Length of product	346.6 mm
Sample volume	50 ml

Product characteristics

Design	with tip, plugged
Label/Print	With print

Color of print/label	violet/black
Graduation	yes

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	30
Piece(s) / case	90
Piece(s) / pallet	6300
Height of inner box	35 mm
Width of inner box	495 mm
Depth of inner box	265 mm
Height of case	149 mm
Width of case	312 mm
Depth of case	409 mm
Volume of case	0.019 cbm
Weight of case	2.46 kg
EAN inner box	4038917394956
EAN case	4038917394949

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT