




## Product specification

### Inoculation loop, 10 µl, PS, blue, sterile



## Product information

Product number	86.1562.010
----------------	-------------

## Product descriptions

Name	Inoculation loop, 10 µl, PS, blue, sterile
Product number	86.1562.010
Product description	Inoculation loop, volume: 10 µl, material: PS, blue, sterile, 10 piece(s)/blister

## Material & colours

Color of product	blue
------------------	------

## Size

Width of product	6 mm
Length of product	200 mm

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

## Packaging

Type of inner packaging	blister
Piece(s) / inner box	500
Piece(s) / case	1000
Piece(s) / pallet	132000
Height of inner box	101 mm
Width of inner box	161 mm
Depth of inner box	216 mm
Height of case	120 mm
Width of case	228 mm
Depth of case	338 mm
Volume of case	0.0092 cbm
Weight of case	1.79 kg
EAN inner box	4038917130875
EAN case	4038917023337

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT