




Product specification

Inoculation loop, 1 µl, PS, white, sterile



Product information

| | |
|----------------|-------------|
| Product number | 86.1567.010 |
|----------------|-------------|

Product descriptions

| | |
|---------------------|---|
| Name | Inoculation loop, 1 µl, PS, white, sterile |
| Product number | 86.1567.010 |
| Product description | Inoculation loop, volume: 1 µl, material: PS, white, sterile, 10 piece(s)/blister |

Material & colours

| | |
|------------------|------------------|
| Product material | Polystyrene (PS) |
| Color of product | white |

Size

| | |
|-------------------|--------|
| Width of product | 4 mm |
| Length of product | 198 mm |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Purity standard | sterile |
| Sterilisation | Gamma irradiation |
| Batched | yes |

| | |
|-------------------------|---------------|
| Minimum order quantity | 1000 |
| Type of inner packaging | blister |
| Piece(s) / inner box | 500 |
| Piece(s) / case | 1000 |
| Piece(s) / pallet | 132000 |
| Height of inner box | 101 mm |
| Width of inner box | 161 mm |
| Depth of inner box | 216 mm |
| Height of case | 120 mm |
| Width of case | 228 mm |
| Depth of case | 338 mm |
| Volume of case | 0.0092 cbm |
| Weight of case | 1.78 kg |
| EAN inner box | 4038917130936 |
| EAN case | 4038917023504 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT