




# Product specification

## Inoculation spreader, PS, blue, sterile



### Product information

Product number	86.1569.001
----------------	-------------

### Product descriptions

Name	Inoculation spreader, PS, blue, sterile
Product number	86.1569.001
Product description	Inoculation spreader, material: PS, blue, for spreading a bacterial suspension, sterile, 1 piece(s)/bag

### Material & colours

Product material	Polystyrene (PS)
Color of product	blue

### Size

Width of product	40 mm
Length of product	185 mm

### Product characteristics

Example of use	for spreading a bacterial suspension
----------------	--------------------------------------

### Purity & certification

Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Gamma irradiation

Batched

yes

## Packaging

Minimum order quantity	500
Type of inner packaging	Bag
Piece(s) / case	500
Piece(s) / pallet	32000
Height of case	215 mm
Width of case	298 mm
Depth of case	398 mm
Volume of case	0.0255 cbm
Weight of case	3.03 kg
EAN inner box	4038917130981
EAN case	4038917023559

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT