



Product specification

Serological pipette, with tip, plugged, 25 ml, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic,



Product information

Product number	86.1685.001
----------------	-------------

Product descriptions

Name	Serological pipette, with tip, plugged, 25 ml, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister
Product number	86.1685.001
Product description	Serological pipette, with tip, plugged, 25 ml, coding colour: red, material: PS, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

Size

Diameter	15 mm
Length of product	345 mm

Product characteristics

Design	with tip, plugged
Color of print/label	red

Graduation	yes
------------	-----

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	200
Piece(s) / inner box	25
Piece(s) / case	200
Piece(s) / pallet	10800
Height of case	290 mm
Width of case	250 mm
Depth of case	420 mm
Volume of case	0.0303 cbm
Weight of case	3.98 kg
EAN inner box	4038917133012
EAN case	4038917024655

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT