



Product specification

Serological pipette, with tip, plugged, 25 ml, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic,



Product information

Product number	86.1685.001
----------------	-------------

Product descriptions

Name	Serological pipette, with tip, plugged, 25 ml, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister
Product number	86.1685.001
Product description	Serological pipette, with tip, plugged, 25 ml, coding colour: red, material: PS, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	15 mm
Length of product	345 mm
Sample volume	25 ml

Product characteristics

Design	with tip, plugged
Label/Print	with printed graduations

Color of print/label	red/black
Graduation	yes

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	25
Piece(s) / case	200
Piece(s) / pallet	10800
Height of case	291 mm
Width of case	250 mm
Depth of case	420 mm
Volume of case	0.0306 cbm
Weight of case	3.98 kg
EAN inner box	4038917133012
EAN case	4038917024655

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT