




Product specification

Demeter pipette, with tip, plugged, 1.1 ml, sterile, 25 piece(s)/bag



Product information

Product number	86.1686.025
----------------	-------------

Product descriptions

Name	Demeter pipette, with tip, plugged, 1.1 ml, sterile, 25 piece(s)/Bag
Product number	86.1686.025
Product description	Demeter pipette, with tip, plugged, 1.1 ml, material: PS, sterile, 25 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	4 mm
Length of product	270 mm
Sample volume	1.1 ml

Product characteristics

Design	with tip, plugged
Label/Print	with printed graduations
Color of print/label	black
Graduation	yes

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Piece(s) / inner box	25
Piece(s) / case	1000
Piece(s) / pallet	108000
Height of case	150 mm
Width of case	250 mm
Depth of case	413 mm
Volume of case	0.0155 cbm
Weight of case	3.94 kg
EAN inner box	4038917133043
EAN case	4038917024662

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT