

Product specification

Power cable 5f MS 1000



Product information

Product number	87.7708.825
----------------	-------------

Product descriptions

Name	Power cable 5f MS 1000
Product number	87.7708.825
Product description	Mains lead with ON/OFF switch for 5 MS 1000

Purity & certification

Product category	no medical device IVD
CE certified	CE
Batched	no

Packaging

Piece(s) / pallet	432
Height of case	100 mm
Width of case	188 mm
Depth of case	188 mm
Volume of case	0.0035 cbm
Weight of case	0.46 kg
EAN inner box	4038917341967
EAN case	4038917341950

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

